


Compiling the SDK in C++ (Using Visual Studio Express 2013)

Installed on local machine:


- Windows SDK 7.1
- Datastream SDK 1.3 (x86) and Datastream SDK 1.3 (x64)
- Visual Studio Express 2013 for Desktop

For Datastream SDK 1.3 (x64)


1. Open Visual Studio Express 2013.
2. Click File -> New Project.
3. Select Win 32 Console Application, and specify a name (currently "Console Application2"). Click OK.




4. Click Next on the following screen.




5. Uncheck Precompiled header and Security Development Lifecycle (SDL) checks. Select Empty Project. Click Finish.



6. Open the folder location for the DataStreamSDK
(C:\Program Files\Vicon\DataStream SDK\Win64\CPP)
7. Open the project folder where this latest project is (default location, for example:
C:\Users\FTsui\Documents\Visual Studio
2013\Projects\ConsoleApplication1\ConsoleApplication1)
8. Copy over the following four files (Client, ViconDataStreamSDK_CPP.dll, ViconDataStreamSDK_CPP and ViconDataStreamSDK_CPPTTest) from the DataStreamSDK folder into the project folder.




9. In VSE 2013, right-click on Header Files and select Add->Existing Item. Select Client.




10. Right-click on Resource Files and select Add->Existing Item. Select ViconDataStreamSDK_CPP.dll and ViconDataStreamSDK_CPP.


11. Right-click on Source Files and select Add->Existing Item. Select ViconDataStreamSDK_CPPTTest. Your Solution Explorer should now look like the following.




12. At the top, change Debug to Release. Change Win32 to Configuration Manager.



13. Change the Active Solution Platform from Win32 to New. Select x64 as the platform. Select OK. Select Close.



14. Change Win32 to x64 so the top should now look like this:




15. Under Build, select Build Solution. The last line of the output window should be:

```
1> Finished generating code
1> ConsoleApplication1.vcxproj -> c:\users\ftsui\documents\visual studio 2013\Projects\ConsoleApplication1\x64\Release\ConsoleApplication1.exe
----- Build: 1 succeeded, 0 failed, 0 up-to-date, 0 skipped -----
```


16. Navigate to the folder where the build is. In this instance, it was in:

"C:\Users\FTsui\Documents\Visual Studio 2013\Projects\ConsoleApplication1\x64\Release".


You will know due to the presence of the Console Application of Type Application. From the DataStreamSDK folder, copy the boost_python-vc90-mt-1_48.dll and boost_thread-vc90-mt-1_48.dll into the folder specified above. Your application folder should look like below now.



17. Now click Local Windows Debugger to start compiling.



Alternatively, you can run it from the folder (above), but you will need to copy in the DLLs like before.



18. The command prompt will start streaming data. To pause, hit the pause key on keyboard. To resume, hit Enter.

```

C:\Users\FTsu\Documents\Visual Studio 2013\Projects\ConsoleApplication1\x64\Release\ConsoleAp...
Children (1):
LeftHumerus
  Static Translation: (0, 0, 0)
  Static Rotation Helical: (0, 0, 0)
  Static Rotation Matrix: (1, 0, 0, 0, 1, 0, 0, 0, 1)
  Static Rotation Quaternion: (0, 0, 0, 1)
  Static Rotation EulerXYZ: (-0, 0, -0)
  Global Translation: (0, 0, 0) True
  Global Rotation Helical: (0, 0, 0) True
  Global Rotation Matrix: (0, 0, 0, 0, 0, 0, 0, 0, 0) True
  Global Rotation Quaternion: (0, 0, 0, 0) True
  Global Rotation EulerXYZ: (0, 0, 0) True
  Local Translation: (0, 0, 0) True
  Local Rotation Helical: (0, 0, 0) True
  Local Rotation Matrix: (0, 0, 0, 0, 0, 0, 0, 0, 0) True
  Local Rotation Quaternion: (0, 0, 0, 0) True
  Local Rotation EulerXYZ: (0, 0, 0) True
Segment #4
Name: LeftHumerus
Parent: LeftClavicle
Children (1):
LeftRadius
  Static Translation: (0, 222.033, 0)
  Static Rotation Helical: (0, 0, 0)
  Static Rotation Matrix: (1, 0, 0, 0, 1, 0, 0, 0, 1)

```

For Datastream SDK 1.3 (x86)

1. It will be the same steps as above except that you will get all the necessary files from:

C:\Program Files (x86)\Vicon\DataStream SDK\Win32\CPP

2. You will not need to change the compiler to x64 in Steps 12-14.
3. Your build will not be under an x64 folder as it is in Step 16. Simply look for the Application File type, and copy in the DLLs (boost_python-vc90-mt-1_48.dll and boost_thread-vc90-mt-1_48.dll). Again, if you want to run the executable from this folder, also copy in ViconDataStreamSDK_CPP.dll and ViconDataStreamSDK_CPP.